

# Best Practices for Sampling and Testing Blood for Lead

Lead is a neurotoxin that has been shown to negatively impact neurological development in children upon exposure at any level. As efforts to reduce blood lead levels in children continue, care must be taken to accurately measure lower concentrations of lead in blood samples. This raises a variety of challenges for specimen collection and laboratory testing practices. High-quality blood lead measurements are necessary to ensure proper case management is provided to individual children and that surveillance leads to proper public health interventions for exposed populations. This fact sheet serves to outline quality improvement guidance for the sampling and measurement of blood lead.

Lead is found in the ambient environment and care must be taken to avoid contamination of specimen collection and testing supplies and to use practices that do not contaminate the specimen throughout the sampling and testing processes.

Potential sources of contamination include:

- The environment in which blood collection is performed.
- Contamination on the skin of the person giving blood (esp. finger-stick)
- Collection supplies which come into contact with the blood sample (e.g., tubes, needles, pipette tips, vials).
- The testing environment.

## **BLOOD SAMPLE COLLECTION: A CRITICAL FIRST STEP**

Minimize contamination of specimens by:

Providing detailed written specimen collection instructions to clinical providers who will be collecting the specimens. These should also include details about specimen handling and transport to ensure the specimens are of appropriate quality for testing.

Using venipuncture samples as preferred specimen type to reduce contamination from the patient's skin.

Using supplies designated or pre-screened for trace-metals and encourage clinical partners to do the same

Screening each lot of sampling and testing supplies for lead contamination, if possible.

Collecting blood samples in a clean environment to prevent contamination, away from known sources of lead (e.g., soil, paint, dust)

Collecting a metals analysis evacuated blood tube first, if collecting multiple tubes, to prevent lead contamination carrying over from tubes not pre-screened for lead

Including field blanks, if collection is part of a biomonitoring study.

Develop specimen acceptance and rejection criteria:

Clotted specimens must be rejected as lead will be preferential found in the red blood cells in the clot

Specimens were collected in the wrong type of tubes.

## THE LABORATORY ENVIRONMENT: KEEP IT CLEAN

Minimize contamination from the ambient environment.

Keep corrugated cardboard and brown recycled paper out of the laboratory.

Use shoe cleaners and tacky mats at entrances to collect dust and dirt.

Provide shoe covers at the lab entrance.

Ensure laboratory processes and practices do not contribute to contamination by:

Training laboratory scientists to use laboratory practices that minimize the likelihood of specimen contamination.

Installing an auto sampler enclosure to prevent airborne contamination of the specimens.

Providing airflow either through filtered (e.g., HEPA or ULPA) fans or by attaching building exhaust.

Testing environmental samples and clinical specimens in different laboratory spaces.

Developing detailed cleaning and decontamination procedures for the laboratory.

## SAMPLE PREPARATION: KEEP IT SAFE, CLEAN, ACCURATE AND PRECISE

Use universal precautions when working with clinical specimens to prevent exposure to blood borne pathogens.

Work within a biological safety cabinet (BSC) to protect both the worker as well as the sample because air entering the cabinet is HEPA filtered before being directed to the work surface.

Prepare analytical standards and reagents in a chemical fume hood.

Wear appropriate personal protective equipment (PPE) including: lab coat, safety glasses, and nitrile gloves.

Use high quality reagents (e.g., double distilled), that have lower background levels of lead. Access to high quality lab water ( $\geq 18.2 \text{ M}\Omega\text{-cm}$ ) is critical.

Acid wash bottles and flasks used for reagents before each use; dedicate to lead testing.

Lot screen tubes and auto sample vials used for testing lead.

Mix blood prior to removing an aliquot using a vortex, a rocker, or pipette pump to ensure homogenous sample.

Inspect blood samples for clots and reject those with clots.

Use a dual-syringe dilutor for consistent and complete volumetric transfers when diluting blood for analysis. Diluent (syringe 1) flushes blood (syringe 2) out of the pipette tip making for a more complete transfer. Initial cleaning (flushing diluent through all syringes and tubing) is important prior to sample preparation.

Check lead background levels of reagents before preparing samples, including diluent before and after the dual-syringe diluter and water that will be used in sample preparation.

## POST-MEASUREMENT QUALITY CONTROL

Analyst and supervisor (or designee): Review all analytical data to verify quality control compliance.

Review transcription for errors if data capture is not used.

Follow your laboratory SOP for confirmatory repeat and or reflex testing for all specimens which initially result in blood lead levels greater than or equal to  $5 \mu\text{g/dL}$ .

## RESOURCES

1. APHL. (October 3, 2019). The importance of lot testing for pediatric lead testing [Webinar]. Retrieved from: <https://aphl.weblinkconnect.com/events/The-Importance-of-Lot-Testing-for-Pediatric-Lead-Testing-00063419-863/details>
2. APHL. (2019). Blood lead testing in public health laboratories [Fact Sheet]. Retrieved from: <https://www.aphl.org/aboutAPHL/publications/Documents/EH-2019Oct-Blood-Lead-Testing-Survey.pdf>
3. Centers for Disease Control and Prevention (CDC). (April 2008). LAMP: CDC's Lead and Multi-element Proficiency Program [Fact Sheet]. Retrieved from: [https://www.cdc.gov/labstandards/pdf/lamp/lamp\\_factsheet.pdf](https://www.cdc.gov/labstandards/pdf/lamp/lamp_factsheet.pdf)
4. CDC. (September 20, 2018). Mission Unleaded: How to test children for lead with maximum accuracy [Video]. Retrieved from: <https://youtu.be/g2p2qREch9g>
5. CDC. (March 2018). Improving the collection and management of human samples used for measuring environmental chemicals and nutrition indicators [Federal Guidance]. Retrieved from: [https://www.cdc.gov/biomonitoring/pdf/Human\\_Sample\\_Collection-508.pdf](https://www.cdc.gov/biomonitoring/pdf/Human_Sample_Collection-508.pdf)
6. Tevis, Denise S., et al. (2018). Assessing the stability of Cd, Mn, Pb, Se, and total Hg in whole human blood by ICP-DRC-MS as a function of temperature and time [Scientific Article]. Clinica Chimica Acta, 485, pp 1-6. Retrieved from: <https://www.sciencedirect.com/science/article/abs/pii/S0009898118302614?via%3Dihub>

This publication was supported by Grant or Cooperative Agreement number #5 NU600E000103, funded by the Centers for Disease Control and Prevention. Its contents are solely the responsibility of the authors and do not necessarily represent the official views of the Centers for Disease Control and Prevention or the Department of Health and Human Services.

© Copyright 2020, Association of Public Health Laboratories. All Rights Reserved.



### ASSOCIATION OF PUBLIC HEALTH LABORATORIES

The [Association of Public Health Laboratories \(APHL\)](http://www.aphl.org) works to strengthen laboratory systems serving the public's health in the US and globally. APHL's member laboratories protect the public's health by monitoring and detecting infectious and foodborne diseases, environmental contaminants, terrorist agents, genetic disorders in newborns and other diverse health threats.

8515 Georgia Avenue, Suite 700, Silver Spring, MD 20910 | 240.485.2745 | [www.aphl.org](http://www.aphl.org)