

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 1 of 21**

1. **PURPOSE:** This procedure describes a standardized laboratory protocol for DNA library preparation of enteric bacterial organisms using the Illumina® DNA Prep kit for subsequent sequencing on an Illumina platform, thus ensuring inter-laboratory comparability of sequencing results and optimizing sequencing data output and quality.

2. **SCOPE:** For use by PulseNet WGS certified laboratorians when preparing libraries from DNA from enteric organisms using the Illumina® DNA Prep kit for sequencing on Illumina platforms for submission of sequencing data to PulseNet. Laboratories may amend this procedure as necessary for use within their laboratories after validation per their laboratory's guidelines.

**NOTE1:** *Instructions for combining Norovirus and Cyclospora amplicons in the same run with the PulseNet organisms are outlined in the appendices PNL35-1 and PNL35-2, respectively.*

**NOTE2:** *This SOP deviates from the protocol provided by Illumina because the original Illumina protocol was optimized for 2 x 150 (300c) sequencing. PulseNet's protocol has been optimized to allow for 2 x 250 (500c) sequencing. If your laboratory follows the original Illumina SOP, the workbooks provided with the PulseNet SOP will not generate accurate calculations.*

3. **DEFINITIONS:**

- 3.1. **BaseSpace:** Illumina cloud-based computing environment for next generation sequencing data analysis, management, and storage, including data sharing
- 3.2. **BLT:** Bead-Linked Transposome
- 3.3. **CD Index:** Combinatorial Dual Index
- 3.4. **CSV:** Comma-Separated Values (file) or Comma Delimited (file)
- 3.5. **DNA:** Deoxyribonucleic Acid
- 3.6. **dsDNA:** Double-Stranded DNA
- 3.7. **EPM:** Enhanced PCR Mix
- 3.8. **Fastq:** A text-based file format for storing both a sequence and its corresponding quality scores
- 3.9. **GHS:** Globally Harmonized System
- 3.10. **HS:** High Sensitivity
- 3.11. **IPB:** Illumina Purification Beads
- 3.12. **LRM:** Local Run Manager
- 3.13. **Mb:** Megabase
- 3.14. **Ng:** Nanogram
- 3.15. **nM:** Nanomolar
- 3.16. **PCR:** Polymerase Chain Reaction
- 3.17. **PHL:** Public Health Laboratory
- 3.18. **PN:** PulseNet
- 3.19. **PPE:** Personal Protective Equipment
- 3.20. **PulseNet Central:** PulseNet team at CDC comprising of the PulseNet Response and Outbreak Management Team ([PulseNet@cdc.gov](mailto:PulseNet@cdc.gov)) and the WGS Core Laboratory Activity ([PulseNetNGSlab@cdc.gov](mailto:PulseNetNGSlab@cdc.gov))
- 3.21. **PulseNet/OutbreakNet SharePoint:** A closed, web-based collaboration application used for communication among PulseNet participants.
- 3.22. **QC:** Quality Control
- 3.23. **RSB:** Resuspension Buffer
- 3.24. **SDS:** Safety Data Sheet
- 3.25. **SOP:** Standard Operating Procedure
- 3.26. **TB1:** Tagmentation Buffer 1
- 3.27. **TSB:** Tagment Stop Buffer
- 3.28. **TWB:** Tagment Wash Buffer

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 2 of 21**

3.29. **UD Index:** Unique Dual Index

3.30. **WGS:** Whole Genome Sequencing

#### 4. RESPONSIBILITIES:

##### 4.1. PulseNet Public Health Laboratory:

4.1.1. Prepare DNA libraries, and QC as necessary, for subsequent WGS.

4.1.2. Re-sequence any isolates that do not meet quality thresholds.

4.1.3. Inform PulseNet Central, as necessary, about any complications with laboratory protocols or suspected reagent issues.

##### 4.2. PulseNet Central:

4.2.1. Perform additional sequence quality analysis in order to provide feedback and troubleshooting support for PHLs as necessary.

4.2.2. Notify PN PHL if any sequences submitted do not meet quality thresholds.

4.2.3. Communicate any suspected reagent issues to PHLs as necessary.

4.2.4. Maintain and review SOPs on a regular basis and post on SharePoint.

#### 5. SAFETY:

5.1. **Biosafety Warning:** This document describes handling of DNA and associated products and does not describe best practices for handling of biological infectious material.

5.2. **Chemical Safety Warning:** Take proper precautions and wear appropriate PPE when handling potentially hazardous chemicals. Ensure that chemicals, spent containers, and unused contents are disposed of in accordance with local and governmental safety standards. **See all relevant SDSs for additional information.**

5.2.1. Illumina® DNA Prep Kit:

5.2.1.1. TSB: GHS Category 1 for eye damage/irritant and is harmful to aquatic life.

5.2.1.2. TB1: GHS Category 4 for acute toxicity (dust/mist), Category 2A for eye irritant and Category 1B for reproductive toxicity. Contains N, N-Dimethylformamide.

5.2.1.3. EPM: GHS Category 4 for acute oral toxicity and Category 1 for specific organ toxicity. Contains tetramethylammonium chloride.

5.2.2. Ethanol is flammable (GHS Flammability Category 2).

#### 6. REAGENTS:

6.1. Illumina® DNA Prep (M) Tagmentation Kit (any of the following):

6.1.1. 96 samples, with IPB, Illumina Cat# 20060059

6.1.2. 24 samples, with IPB, Illumina Cat# 20060060

Kit components:

- **Beads + Buffers.** Store at 15 to 30°C. **NOTE:** *the outside label lists storage at 2-8°C but all component vials list storing at 15 to 30°C.*

**IPB (store upright)**

TSB

TWB

- **PCR + Buffers.** Store at -25 to -15°C.

RSB

TB1

EPM

- **Tagmentation (M) Beads.** Store at 2-8°C.

BLT (store upright)

6.2. Indexes: choose either 6.2.1, 6.2.2 or 6.2.3.

**NOTE1:** *The index kit(s) needed depends on the sample throughput of the laboratory and the level of*

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 3 of 21**

*multiplexing desired. It is recommended not to use the same index pair in two consequent runs on the same instrument to prevent carry-over from run to run.*

**NOTE2:** *The CD and UD indexes **cannot** be mixed on the same run.*

**NOTE3:** *The IDT for Illumina DNA/RNA indexes and the Illumina DNA/RNA indexes **cannot** be mixed on the same run.*

**NOTE4:** *The UD indexes have TEN base pair adapters; therefore, the number of index read cycles on the instrument must be adjusted to 10 (from 8 for CD indexes).*

- 6.2.1. Nextera DNA CD Indexes. Store at -25°C to -15°C. 96 Dual Index, plate format (96 samples, Illumina Cat# 20018708)
- 6.2.2. IDT for Illumina DNA/RNA UD Indexes. Store at -25°C to -15°C.

**NOTE:** *These index kits will be discontinued: last order date January 29<sup>th</sup>, 2025 or until inventory is depleted.*

  - 6.2.2.1. Set A, Tagmentation (96 indexes, 96 samples, Illumina Cat# 20027213)
  - 6.2.2.2. Set B, Tagmentation (96 indexes, 96 samples, Illumina Cat# 20027214)
  - 6.2.2.3. Set C, Tagmentation (96 indexes, 96 samples, Illumina Cat # 20042666)
  - 6.2.2.4. Set D, Tagmentation (96 indexes, 96 samples, Illumina Cat # 20042667)
- 6.2.3. Illumina DNA/RNA UD Indexes. Store at -25°C to -15°C.
  - 6.2.3.1. Set A, Tagmentation (96 indexes, 96 samples, Illumina Cat # 20091654)
  - 6.2.3.2. Set B, Tagmentation (96 indexes, 96 samples, Illumina Cat # 20091656)
  - 6.2.3.3. Set C, Tagmentation (96 indexes, 96 samples, Illumina Cat # 20091658)
  - 6.2.3.4. Set D, Tagmentation (96 indexes, 96 samples, Illumina Cat # 20091660)
- 6.3. Ethanol, molecular grade, 95-100% (Fisher Scientific Cat# BP2818-500 or equivalent)
- 6.4. Ethanol, lab-grade, 70% or equivalent for disinfection purposes (Fisher Scientific Cat# 04-355-305 or equivalent)
- 6.5. Water, Molecular grade (Fisher Scientific Cat# BP24701 or equivalent)
- 6.6. Invitrogen™ Qubit™ dsDNA HS Assay Kit: choose either 6.6.1. or 6.6.2.
  - 6.6.1. Concentrated assay reagent kit (100 assays, Fisher Scientific Cat# Q32851 **OR** 500 assays, Fisher Scientific Cat# Q32854) with the following components:
    - dsDNA HS Reagent, Component A (room temperature, protect from light)
    - dsDNA HS Buffer, Component B (room temperature)
    - dsDNA HS Standard #1, Component C ( $\leq 4^{\circ}\text{C}$ )
    - dsDNA HS Standard #2, Component D ( $\leq 4^{\circ}\text{C}$ )
  - 6.6.2. 1x ready-to-use assay reagent kit (100 assays, Fisher Cat# Q33230 **OR** 500 assays Fisher Scientific Cat# Q33231) with the following components (2-8°C):
    - dsDNA HS Working Solution, Component A (room temperature, protect from light)
    - dsDNA HS Standard #1, Component B ( $\leq 4^{\circ}\text{C}$ )
    - dsDNA HS Standard #2, Component C ( $\leq 4^{\circ}\text{C}$ )

## **7. SUPPLIES:**

- 7.1. 96 well PCR plates, skirted, hard shell low profile, thin-wall (BioRad Cat# HSP-9601 or equivalent)
- 7.2. Deepwell storage “MIDI” plates, 96 well (Fisher Scientific Cat# AB-0859 or equivalent) – **Optional**
- 7.3. Ice
- 7.4. Microcentrifuge tubes, 1.5 ml, sterile (Fisher Scientific Cat# 05-408-129 or equivalent)
- 7.5. Microseal B adhesive seal (BioRad Cat# MSB-1001 or equivalent)
- 7.6. Microseal F adhesive foil seal (BioRad Cat# MSF-1001 or equivalent) – **Optional**
- 7.7. Pipette tips, sterile, filtered: 20  $\mu\text{l}$ , 200  $\mu\text{l}$  and 1000  $\mu\text{l}$  volumes (Rainin Cat# 30389225, 30389239 & 30389212 or equivalent)
- 7.8. Qubit™ Assay Tubes (Fisher Scientific Cat# Q32856 or equivalent (clear, thin-wall 0.5-ml PCR tubes)).
- 7.9. Solution basins, sterile (Fisher Scientific Cat# 13-681-504 or equivalent)

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 4 of 21**

**8. EQUIPMENT:**

- 8.1. Ice buckets/containers
- 8.2. Invitrogen™ DynaMag™-96 Side Skirted Magnet (Fisher Scientific Cat# 12-027 or equivalent)
- 8.3. Microcentrifuge for quick spins
- 8.4. Micropipettes, capable of volumes from 1 µl to 1000 µl. Single and multichannel (20 µl and 100µl volumes).

**NOTE:** *Two sets of pipettes are suggested, one for working with pre-amplified product and reagents and another set for working with post-PCR amplified product and reagents.*

- 8.5. Microplate centrifuge
- 8.6. Qubit™ 4.0 Fluorometer (older versions 2.0 or 3.0), or equivalent for quantification of dsDNA
- 8.7. Thermal cycler with heated lid, compatible with 96-well 0.2 ml PCR plate
- 8.8. Vortex

- 9. PROCEDURE:** The Illumina DNA Prep checklist (PNL35.W3) may be used during library preparation in the laboratory as a record of what steps have been completed. This checklist will also calculate master mix volumes. A quick guide (PNL35.JA1) outlining the library preparation method is also available as a laboratory reference for the library preparation procedure. The DNA Prep workbook (PNL35.W1) may be used for planning sequencing runs (on the iSeq, MiniSeq, MiSeq or NextSeq), including index assignment and may also be printed and used in the laboratory to record reagent lots and other run-specific information. The Index tracking workbook (PNL35.W4) may be used for assignment of indexes and for inventory purposes.

**NOTE1:** *The only safe stopping points in the procedure are after the amplification step (9.4.), after the clean-up of the amplified libraries (9.5.) and after dilution to 4 nM (9.6.).*

**NOTE2:** *All multichannel mixing steps in this protocol may be replaced by the use of a plate shaker (1600 rpm for 1 minute) if desired.*

**NOTE3:** *Make sure the vortex is turned down once cells have been lysed during DNA extraction and for all steps during library prep.*

**NOTE4:** *It has been documented and observed that cluster density can be negatively affected by the use of ammonium-based cleaning products near sequencing equipment, including lab benches and pipets used for library preparation. Do not use quaternary ammonium compounds or wipes near or on sequencing equipment!*

**9.1. (Optional) Prepare the “Library Prep” Tab of the Illumina DNA Prep Workbook**

**NOTE:** *Do NOT delete rows or columns in the workbook PNL35.W1! There are many formulas and look-up tables that could be disturbed by doing so. There are numerous rows available for sample information; any rows which are not going to be used may be hidden by using the “Hide” feature in Excel. The workbook is designed with the following color scheme:*

- *White fields should be filled in*
- *Dark gray fields are optional*
- *Blue fields contain formulas, which will auto populate, and should not be altered*
- *Yellow index wells indicate that duplicate indexes were selected*
- *Red index wells indicate that CD indexes were selected with v2 and/or v3 indexes*

9.1.1. Enter the Run ID using the format PulseNet lab ID-Instrument ID-run start date: labID-MXXXX-YYMMDD. (e.g., CDC-M3235-240512).

9.1.2. Enter the Library Prep Date, Library Prep Technician, Sequencing Kit Type/Chemistry from the drop-down menu, Sequencing Date, and Sequencing Technician.

9.1.3. Enter the State Keys (the ID entered in the PulseNet database “Key” field).

**NOTE: Important!** *Fastq file names will be assigned on the Sample Sheet tab according to this State Key and the Run ID. These fields will be concatenated to create a unique prefix*

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

Doc. No. PNL35

Ver. No. 05

Effective Date:

Page 5 of 21

for the resulting fastq files (e.g., Sample1-CDC-M3235-240512)

9.1.4. Enter the Genome Size Estimate (based on Table 1 below) for each organism in Column E of the workbook.

Organism	Estimated Genome Size (Mb)	Minimum Coverage
<i>E. coli &amp; Shigella spp.</i>	5	40x
<i>Salmonella spp.</i>	5	30x
<i>Vibrio cholerae</i>	4	40x
<i>Vibrio spp.</i>	5	40x
<i>Listeria monocytogenes</i>	3	20x
<i>Campylobacter spp.</i>	1.6	20x
<i>Cronobacter spp.</i>	4.5	40x
<i>Yersinia spp.</i>	5	40x

**Table 1.** Estimated genome size (in Mb) by organism

9.1.5. Confirm that the number of isolates on the run is appropriate for the capacity of the instrument and reagent kit to be used. The recommended DNA load for each of the Illumina sequencing platforms and cartridges validated for use by PulseNet can be found in Table 2 below. Individual laboratories with experience and consistent library concentrations and fragment sizes may be able to exceed the recommended DNA load while still reaching the required minimum coverage levels (Table 1) for all isolates on the run. The sum of the genome sizes (in Mb) for the samples on the run will be displayed on the workbook and is equal to the estimated DNA load of the run.

**NOTE1:** *The lower recommended DNA load for runs containing E. coli, Shigella spp., Vibrio spp., Cronobacter spp. or Yersinia spp. is due to the fact that these organisms require a higher coverage (40x) for downstream analyses compared to the other organisms (20x or 30x).*

*Therefore, more sequence reads/data generated in the sequencing run is required for these organisms to meet coverage requirements.*

**NOTE2:** *PulseNet never accepts, under any circumstances, sequences generated with cycling parameters below 300 cycles (2 x 150).*

Sequencing Kit (cycles)	DNA Load (Mb) Runs Without 40x Organisms	DNA Load (Mb) Runs Containing 40x Organisms
<b>MiSeq</b>		
v2 300	90	90
v2 500	100	100
v3 600	200	175
Micro (300)	35	30
Nano (500)	13	13
Nano (300)	13	10
<b>MiniSeq</b>		
Mid Output	60	60
High Output	100	100
<b>iSeq</b>		
iSeq v2	35	25

**Table 2.** Estimated DNA load capacity (in Mb) for Illumina reagent kits.

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 6 of 21**

9.1.6. Determine which indexes will be used and enter the index plate well position (from the index plate). The index plate well may be entered manually, selected from the drop-down option or copied over from the index tracking worksheet.

**NOTE1:** *It is recommended not to use the same index pairs within two consecutive runs on the same sequencer to reduce the amount of carryover. The index tracking worksheets, PNL35.W4 may be used to reserve/keep track of used and available indexes.*

**NOTE2:** *CD and UD indexes cannot be combined on the same run. It is also not recommended to combine the IDT for Illumina DNA/RNA UD indexes with the Illumina DNA/RNA UD indexes.*

**NOTE3:** *it is highly recommended to use the UD indexes for this procedure; they have less index hopping in the instruments utilizing patterned flow cells (NextSeq, MiniSeq, iSeq) and are color balanced because all index sequences are unique.*

9.1.7. Enter the volume of input DNA (recommended 10 µl) to be added to each well at the beginning of the library preparation process.

**NOTE1:** *PulseNet has standardized the starting volume of DNA to 10 µl; however, the recommended quantity of input DNA per Illumina is 100-500 ng. Individual laboratories may adjust input DNA volumes to ensure DNA falls within this range if desired.*

**NOTE2:** *The minimum volume of input DNA is 2 µl. If DNA is too concentrated, perform a dilution to bring input DNA volume above 2 µl and proceed.*

9.1.8. The DNA Prep tab of the workbook may now be printed for use in the laboratory.

**9.2. Dilute and Tagment Input DNA: DNA is fragmented and tagged with the adapter sequences and bound to the BLT during these steps. It is important that beads are well suspended at ALL steps in the procedure.**

**NOTE:** *Ensure that DNA going into library preparation has been assessed for quality. The 260/280 value should be between 1.75 and 2.05. See PNL33 for more information.*

9.2.1. Allow BLT (from the refrigerator) and TB1 (from freezer) to come to room temperature.

**NOTE:** *Ensure that BLT is never frozen and is stored upright at all times so that the beads always remain submerged in the buffer.*

9.2.2. Label a 96-well PCR plate, or equivalent, with Run ID.

9.2.3. Add molecular grade water to each sample well (Total assay input volume is 30 µl; water volume is 20 µl if 10 µl of DNA will be added; volume is automatically calculated and displayed on the Workbook).

9.2.4. Add DNA to the molecular grade water (per volume entered in the Workbook, typically 10 µl) and mix well by gently pipetting approximately 5-10 times.

9.2.5. Vortex BLT for a **minimum** of 10 seconds and verify proper suspension of beads; repeat if necessary. Do not centrifuge.

9.2.6. Vortex TB1 to mix and perform a quick spin.

9.2.7. Prepare tagmentation master mix:

Reagent	Volume per Sample
TB1	10 µl
BLT	10 µl

**Table 3:** Reagent volumes per sample for tagmentation master mix

**NOTE:** *The number of samples being prepared may be entered into the Illumina DNA Prep checklist, PNL35.W3 (“checklist”) and master mix volumes will be auto calculated with surplus of 2 samples already included in the calculation.*

9.2.8. Vortex the tagmentation master mix well.

9.2.9. Add 20µl of tagmentation master mix to each sample well.

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**


**Ver. No. 05**

**Effective Date:**

**Page 7 of 21**

**NOTE:** *If many samples (i.e., greater than 8) are being prepped at once, it may be necessary to re-vortex the tagmentation master mix occasionally to ensure that it does not settle during this step.*

- 9.2.10. Pipet gently to mix well, resuspending the beads (using a multichannel pipette for mixing is recommended). Inspect the wells to ensure that samples and master mix are homogeneous before incubation. See Figure 1 below for additional information.



**Figure 1.** Example of homogeneous and non-homogeneous sample wells. Wells 1, 4, 6, 7 and 8 are well mixed. There should be no color striations (“sunset”) from the bottom to the top of the wells. Wells 2, 3 and 5 display this settling and should be re-mixed prior to incubation.

- 9.2.11. Seal the plate with Microseal B (or equivalent) and incubate the plate at 55°C for 15 minutes, followed by a 10°C hold on a thermal cycler with the lid heated at 100°C (volume is 50 µl).

**NOTE:** *It is recommended to pre-program a thermal cycler for this purpose.*

**9.3. Clean Up Tagmented DNA: DNA (now tagged with adapters and bound to the BLT) will be washed prior to amplification.**

**NOTE:** *A deep-well MIDI plate may be used for the washing steps below if desired. This may improve formation of the pellet and decrease cross-contamination risks. The samples will need to be transferred to a 96-well skirted PCR plate for thermal cycler steps.*

- 9.3.1. Ensure that TSB is at room temperature and check that no precipitate is present (if precipitated, warm at 37°C for up to 10 minutes or until dissolved and vortex).

- 9.3.2. Add 10 µl of TSB to each sample.

- 9.3.3. Pipet gently to mix. Inspect the wells to ensure that samples and master mix are homogeneous prior to incubation. Refer to Figure 1.

- 9.3.4. Seal the plate with Microseal B or equivalent.

- 9.3.5. Incubate at 37°C for 15 minutes and hold at 10°C on a thermal cycler with a lid pre-heated to 100°C.

**NOTE:** *It is recommended to pre-program a thermal cycler for this purpose.*

- 9.3.7. After samples reach 10°C, remove the plate from the thermal cycler, and perform a quick spin.

- 9.3.8. Place on a magnet for 3 minutes (or until beads form a tight pellet).

**NOTE:** *If a tight pellet is not formed after 3 minutes, gently pipet the samples that haven't formed a pellet and allow to sit on the magnet for up to 3 minutes.*

- 9.3.9. Remove and discard supernatant.

- 9.3.10. Remove the plate from the magnet and add 100 µl TWB directly to the pellet.

**NOTE:** *A solution basin and multichannel pipet may be used for TWB washes if desired.*

- 9.3.11. Pipet gently (to avoid foaming) to mix until beads are fully resuspended.

**NOTE:** *Pipet with care. Some foaming may occur, but excessive foaming of TWB should be avoided, as foaming can lead to sample loss.*

- 9.3.12. Place back on the magnet for 3 minutes (or until beads form a tight pellet).

- 9.3.13. Remove and discard the supernatant.

- 9.3.14. Remove from the magnet and gently add 100 µl of TWB directly to the pellet.

- 9.3.15. Pipet gently to mix (avoiding foaming) until beads are fully resuspended.

- 9.3.16. Place back on the magnet for 3 minutes (or until beads form a tight pellet).

- 9.3.17. Remove and discard supernatant.

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 8 of 21**

- 9.3.18. Remove from the magnet, add 100µl of TWB directly to the beads and gently pipet to resuspend.  
9.3.19. Place back on the magnet for 3 minutes, allowing TWB to remain in the wells (to prevent drying of beads) while proceeding to amplification steps.

**9.4. Amplify Tagmented DNA: The indexes are added to the tagmented DNA, and amplification occurs. DNA will be released from the BLT beads during PCR.**

- 9.4.1. Thaw EPM on ice and thaw indexes at room temperature  
9.4.2. Vortex and quick spin the thawed EPM briefly prior to use.  
9.4.3. Prepare PCR master mix:

Reagent	Volume per Sample
EPM	20 µl
Molecular grade water	20 µl

**Table 4:** Reagent volumes per sample for PCR master mix.

**NOTE:** *The number of samples being prepared may be entered into the Illumina DNA Prep checklist and master mix volumes will be auto calculated with surplus of 2 samples already included in the calculation.*

- 9.4.4. Vortex and quick spin the PCR master mix.  
9.4.5. Remove TWB from wells.  
**NOTE:** *Removal of TWB is crucial, as it can impede PCR. Using a small volume pipette to remove all residual TWB is recommended. Any foam remaining on the wells will not negatively impact the library.*  
9.4.6. Remove from the magnet and immediately add 40 µl of PCR master mix to each sample.  
**NOTE:** *It is recommended to ensure that each pellet is submerged or wetted by the master mix prior to moving on to the next well and to add master mix to all wells prior to mixing in order to avoid beads drying out.*  
9.4.7. Gently pipet to mix well, ensuring resuspension of the pellet.  
9.4.8. Add indexes to samples:  
9.4.8.1. Each index plate well contains a unique index pair in a single-use volume. Add 10 µl of appropriate index pair from the index plate well to each sample well as indicated by the workbook (PNL35.W1).  
9.4.8.2. After wells have been used/pierced they **must be sealed** using lab tape or Microseal F or equivalent to prevent index cross-contamination.  
**NOTE:** *It is recommended to pierce the foil of the desired well on the index plate with a new pipette tip, then to use a fresh pipette tip to remove the indexes from the wells.*  
9.4.9. Pipet a minimum of 10 times to mix.  
9.4.10. Inspect the wells to ensure that all samples are homogeneous. This step is crucial for effective indexing and subsequent product recovery. Refer to Figure 1.

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 9 of 21**

9.4.11. Seal the plate with Microseal B or equivalent, place onto thermal cycler, and run the following pre-programmed settings with a heated lid (100°C):

Step 1: 68°C for 3 minutes

Step 2: 98°C for 3 minutes

Step 3: 5 cycles

98°C for 45 seconds

62°C for 30 seconds

68°C for 2 minutes

Step 4: 68°C for 1-minute

Step 5: Hold at 10°C

Total volume: 50 µl

9.4.12. Centrifuge plate at 280 x g for 1 minute.

**NOTE:** *This is a safe stopping point. The plate may be sealed with Microseal B or equivalent and stored at 2-8°C or left on the thermal cycler at 10°C for up to 3 days.*

**9.5. Clean up and Size Select Amplified Libraries: This dual bead clean-up procedure is for purification and size selection of libraries. Target fragment sizes are 800-1000 bp.**

**NOTE1:** *The steps listed below are critical for efficient size selection, product recovery, and thus cluster generation and sequencing. Always check pipette tips for correct volumes and ensure that no beads have accidentally been aspirated on steps where supernatant is being removed. If beads have accidentally been aspirated or the bead pellet has been disturbed, allow the pellet to re-form and repeat the step. It is also important to ensure that beads are well suspended at all “pipet to mix steps” as described below.*

**NOTE2:** *The target fragment size of 800-1000 bp has been optimized for the PulseNet protocol (instead of the 500 bp fragment size targeted by the original Illumina DNA Prep procedure) in order to better facilitate the use of 500 cycle sequencing chemistry.*

9.5.1. Bring RSB to room temperature (from frozen) and vortex to mix.

**NOTE:** *RSB is in a 50 ml conical. In order to facilitate faster thawing, it is recommended to aliquot and freeze back smaller volumes after opening for the first time*

9.5.2. If the plate was retrieved from cold storage, centrifuge plate at 280 x g for 1 minute.

9.5.3. Place plate on the magnet for 5 minutes (or until beads have formed a tight pellet).

**NOTE:** *If any wells do not have tight pellets after 5 minutes, gently resuspend them and allow time for a tight pellet to form.*

9.5.4. Transfer 45 µl of supernatant (now containing the DNA) to new wells.

**NOTE1:** *A MIDI plate may be used for steps 9.5.5. through 9.5.26. if desired. This may improve pelleting of the beads and reduce the tendency of the beads to be pulled off of the magnet during removal of the final product during final elution step.*

**NOTE2:** *If sample loss has occurred during tagmentation and PCR, it is important to maintain the appropriate ratio of IPB master mix to DNA (1.88x) for recovery of appropriate library lengths. Measure and transfer as much of the resulting supernatant as possible and add the corresponding volume of IPB master mix from the Table 5 below. Do not adjust subsequent volume of stock IPB later in the protocol, as this will negatively affect product recovery yield. If less than 30 µl of sample is available for transfer, do not proceed with size selection but repeat the library prep instead.*

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 10 of 21**

Volume of DNA recovered after tagmentation (µl)	Volume of IPB Master mix to be added (µl)
40	75.2
35	65.8
30	56.4

**Table 5.** Adjusted IPB master mix volumes in the case of sample loss.

9.5.5. Vortex and invert IPB stock several times to fully resuspend the beads.

9.5.6. Prepare **IPB master mix**:

Reagent	Volume per Sample
IPB	40.8 µl
Molecular grade water	44.2 µl

**Table 6:** Reagent volumes per sample for IPB master mix

**NOTE:** *The number of samples being prepared may be entered into the DNA Prep checklist and master mix volumes will be auto calculated with surplus of 2 samples already included in the calculation.*

9.5.7. Remove sample plate from the magnet.

9.5.8. Vortex IPB master mix thoroughly and add 85µl to each sample (or appropriate volume from Table 5 above, if applicable).

9.5.9. Pipet to mix a minimum of 10 times; using a multichannel pipet is recommended.

**NOTE:** *Use caution when mixing as the volume will be > 100 µl.*

9.5.10. Incubate at room temperature for 5 minutes.

9.5.11. Place on the magnet for 3 minutes (or until beads form a tight pellet).

9.5.12. During incubation, vortex the stock IPB.

9.5.13. After incubation, with the plate still on the magnet, transfer 125 µl of supernatant (containing DNA) to new wells (longest DNA fragments are being left behind).

**NOTE:** *If volumes of IPB were changed due to product loss, the volume of supernatant transferred to new wells will be lower. Transfer as much of the supernatant as possible.*

9.5.14. Remove the plate from the magnet and add 15 µl of well-suspended **stock IPB** to the supernatant (undiluted beads to DNA ratio 0.12).

**NOTE:** *Be sure to add from the stock tube of undiluted IPB, NOT the previously made IPB master mix. Accidentally using the IPB master mix will dramatically decrease yield.*

9.5.15. Gently pipet a minimum of 10 times to mix.

9.5.16. Incubate at room temperature for 5 minutes.

9.5.17. Dilute enough fresh 80% ethanol for all samples:

Reagent	Volume per sample	Example: 20 samples
100% ethanol	0.4 ml	8 ml
Molecular grade water	0.1 ml	2 ml

**Table 7:** Reagent volumes per sample for 80% ethanol.

9.5.18. Place sample plate on magnet for 3 minutes (or until beads form a tight pellet and supernatant clears).

9.5.19. Remove and discard supernatant (desired libraries are bound to the beads; shortest libraries in the supernatant are being discarded).

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 11 of 21**

9.5.20. Perform the steps below (9.5.20.1. through 9.5.20.3.) twice, for a total of two washes:

9.5.20.1. While the plate is on the magnet, add 170  $\mu$ l of fresh 80% ethanol.

**NOTE1:** *Do not add directly to the bead pellet, do not mix, and do not remove from the magnet during the wash steps.*

**NOTE2:** *A solution basin may be used for the ethanol if desired.*

9.5.20.2. Incubate for 30 seconds.

9.5.20.3. Remove and discard supernatant.

9.5.21. After the second wash, use a small volume pipette to remove excess ethanol, if necessary.

9.5.22. Allow beads to air dry for 3-5 minutes.

**NOTE:** *Do not allow beads to over-dry or crack. It is advisable to err on the shorter side of drying time. If cracking is observed, immediately resuspend beads as described below regardless of drying time.*

9.5.23. Remove from magnet and add 32  $\mu$ l of RSB. It is recommended to add RSB quickly to each sample pellet, then go back and mix with a multichannel pipet after all pellets have been wetted.

9.5.24. Pipet gently and thoroughly to mix.

9.5.25. Incubate at room temperature for 2-5 minutes.

**NOTE:** *Longer incubation (5-10 minutes) is preferred for optimal yield and recovery of longer libraries. If yield is low (less than 10 ng/ $\mu$ l), RSB incubation time may be increased to 8 minutes.*

9.5.26. Place on magnet for 3 minutes (or until supernatant is clear).

9.5.27. Transfer 30  $\mu$ l of the supernatant into new wells (or into wells on a new plate) – this is the final product.

**NOTE1:** *If a skirted PCR plate is used for this step, the pellet may slide away from the magnet. To prevent aspirating beads, decreasing the recovery volume to 25  $\mu$ l is acceptable.*

**NOTE2:** *A quantification step (Qubit dsHS kit or equivalent) may be performed at this point on the individual libraries and results may be recorded in the Workbook if desired. See PNL33 for specific instructions on operating the Qubit. Checking post-clean up library concentrations is recommended for new users of the Illumina DNA Prep kit and for troubleshooting purposes. Ideal library yield is 10-20 ng/ $\mu$ l. Sequencing libraries successfully with yields below 10 ng/ $\mu$ l is possible but indicates inefficient library preparation. Library concentrations significantly higher than 20 ng/ $\mu$ l also indicate a likely error in the library prep process and should be repeated.*

**NOTE3:** *The quality of the individual libraries may also be checked using fragment analysis at this point, particularly for new users of the Illumina DNA Prep kit, for troubleshooting purposes or when multiplexing libraries from multiple preparation batches for a high throughput sequencing instrument. The average fragment size should be 800-1000 bp and the fragment distribution graph should display a tight well-formed peak. Libraries with very small or large average fragment sizes, very broad peaks or double peaks should be re-prepped.*

**NOTE4:** *This is a safe stopping point. The plate may be sealed with Microseal B or equivalent and stored at -20°C for up to 30 days per Illumina recommendation. Longer term storage is possible but repeated freeze-thaw cycles of the libraries will result in loss of yield.*

## 9.6. Pool Libraries

9.6.1. Pool 5  $\mu$ l of each library into a new well, with the exceptions noted below in Table 8, and pipet well to mix.

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 12 of 21**

Organism	Instrument	Pooling volume (µl)
<i>Campylobacter</i>	MiSeq, MiniSeq	2.5
<i>Campylobacter</i>	iSeq	2.0
<i>Listeria</i>	iSeq	2.0

**Table 8.** Exceptions to the 5 µl pooling volume for the Illumina DNA Prep libraries

**NOTE:** Lower pooling volumes are required to achieve a more balanced run. This will help prevent over-coverage due to small genome size and lower coverage requirements.

9.6.2. Quantify the pool, using the Qubit dsHS kit or equivalent. See SOP PNL33 for instructions on operating the Qubit.

9.6.3. Optional: perform fragment analysis to determine the average fragment size of the pool.

9.6.4. Enter the concentration (ng/µl) for the pool in the appropriate cell of the Workbook.

9.6.5. Calculate the molarity (nM) of the pool:

$$(\text{Qubit reading ng/}\mu\text{l} / (660\text{g/mol} \times 1000\text{bp})) \times 10^6$$

**NOTE1:** The workbook will automatically calculate and display this value.

**NOTE2:** The formula in the workbook is based on an assumed average fragment size of 1000 base pairs (in blue above) which is optimal for 500 cycle sequencing (and deviates from the original Illumina SOP). If laboratories generate fragment sizes after library prep which vary greatly from 1000 bp in length, they may adjust this formula accordingly, as this will make the determination of the library concentration and dilution to loading concentration more accurate. For example, if fragment analysis reveals average library lengths at around 800 bp, the formula would be the following:  $(\text{Qubit reading ng/}\mu\text{l} / (660\text{g/mol} \times 800\text{bp})) \times 10^6$  and therefore the molarity calculation cell would need to be changed to “ $= (\text{PoolConcentration} / (660 * 800)) * 10^6$ ”. Contact [PulseNetNGSLab@cdc.gov](mailto:PulseNetNGSLab@cdc.gov) for assistance if necessary.

9.6.6. Calculate the volume (µl) of pool necessary to generate 50 µl of a 4 nM pool:

$$(200 / \text{Molarity of pool}) = \text{volume of pool for dilution}$$

**NOTE:** The workbook will automatically calculate and display this value.

9.6.7. Calculate the volume (µl) of RSB diluent required:

$$50 - \text{volume of pool} = \text{volume of RSB required}$$

**NOTE:** The workbook will automatically calculate and display this value.

9.6.8. In a new well, dilute the pool to 4 nM by adding the calculated volume of library pool to the calculated volume of RSB.

**NOTE1:** The pooled libraries are now ready for sequencing, or the plate may be sealed with Microseal B or equivalent and stored at -20°C for no more than 48 h. If continuing with sequencing, proceed to the appropriate instrument sequencing SOP for instructions (PNL38, PNL39, PNL40) to denature and/or dilute the pool to the proper loading concentration. The SOP for the NextSeq sequencing is still under development but the appendix PNL35-3 outlines the steps needed to utilize the sample import template tab in the PNL35.W1 to generate a sample sheet for the NextSeq.

**NOTE2:** To include Norovirus or Cyclospora to the PulseNet sequencing run, refer to Appendices PNL35-1 and PNL35-2, respectively, for instructions on how to spike in Norovirus and Cyclospora libraries to the 4 nM bacterial library pool.

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

Doc. No. PNL35


Ver. No. 05

Effective Date:


Page 13 of 21

**10. FLOW CHARTS:**


**10.1. PulseNet WGS Workflows:**



**10.2. Illumina DNA Prep Library Preparation Workflow:**



**10.3. Detailed Visualization of the Illumina DNA Prep Library Preparation Workflow**




**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 14 of 21**



**11. RELATED DOCUMENTS:**

Document Number	Title
PNL33	DNA Extraction and QC SOP
PNL38	Sequencing on the MiSeq SOP
PNL39	Sequencing on the MiniSeq SOP
PNL40	Sequencing on the iSeq SOP
PNQ07	Illumina Sequence Data QC SOP
PNL35.W1	Illumina DNA Prep Workbook, 96 CD/UD Indexes, separate sample sheets for MiSeq LRM, MiniSeq LRM and iSeq LRM, and a Sample Import Template for NextSeq LRM
PNL35.W3	Illumina DNA Prep Checklist Workbook
PNL35.W4	Illumina DNA Prep Index Tracking Template Worksheet
PNL35.JA1	Illumina DNA Prep Quick Guide Job Aid
VGB.NV.METHOD.002	RT-PCR Protocol for Amplification of Full-Length Norovirus GII Genomes
AMD.DR.C.001	State Cyclospora Outbreak Genotyping SOP

**12. REFERENCES:**

- 12.1. Illumina, Inc. Index Adapters Pooling Guide. (Doc.# 1000000041074 v13). February 2024. <https://support-docs.illumina.com/SHARE/IndexAdaptersPooling/Content/SHARE/FrontPages/IndexAdapterPooling.htm>
- 12.2. Illumina, Inc. Illumina DNA Prep Checklist (Doc.# 1000000033561 v05). June 2020. [https://support.illumina.com/content/dam/illumina-support/documents/documentation/chemistry\\_documentation/illumina\\_prep/illumina-prep-checklist-1000000033561-05.pdf](https://support.illumina.com/content/dam/illumina-support/documents/documentation/chemistry_documentation/illumina_prep/illumina-prep-checklist-1000000033561-05.pdf)
- 12.3. Illumina, Inc. Illumina DNA Prep Reference Guide (Doc.# 1000000025416 v09). June 2020. [https://support.illumina.com/content/dam/illumina-support/documents/documentation/chemistry\\_documentation/illumina\\_prep/illumina-dna-prep-reference-guide-1000000025416-09.pdf](https://support.illumina.com/content/dam/illumina-support/documents/documentation/chemistry_documentation/illumina_prep/illumina-dna-prep-reference-guide-1000000025416-09.pdf)
- 12.4. Graphic 1 (basic design): Ko, Julie. PCR Tube. The Noun Project.

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 15 of 21**

<https://thenounproject.com/search/?q=pcr&i=1702287>

12.5. Illumina, Inc. Customer Notification. November 10, 2021.

<https://support.illumina.com/bulletins/2019/09/impact-of-ammonium-based-cleaning-products-on-sequencing-run-per.html>

### 13. CONTACTS:

13.1. PulseNet NGS Lab troubleshooting account: [PulseNetNGSLab@cdc.gov](mailto:PulseNetNGSLab@cdc.gov)

### 14. AMENDMENTS:

14.1. **01/28/19:** New Document

14.2. **04/09/20:** Added notes throughout based on feedback and PHL experience, added Table 7, added workflows and iSeq/MiniSeq information. Made several appendices individual workbooks & added quick guide.

14.3. **03/30/2021:** Updated the SOP throughout with the new name for the Nextera DNA Flex kit: Illumina® DNA Prep. Updated catalog numbers and added the IDT for Illumina UD indexes. Updated the references section. Added appendix PNL35-1 for instructions to make combined runs with PulseNet organisms and viral amplicons and added a workflow chart for the procedure.

14.4. **02/01/2023:** Updated catalog numbers and added notes for library prep and index kits to be obsoleted by Illumina. Added language explaining all significant deviations from the Illumina SOP and language facilitating deviations from this SOP, such as using a higher DNA load capacity than the minimum recommended by PulseNet. Added library pooling instructions for iSeq. Added PNL35-2 for instructions to make combined runs with PulseNet organisms and *Cyclospora* amplicons and added a workflow chart for the procedure. Updated the workbooks to reflect the changes in the index kits (PNL35.W1, PNL35.W4), added the MiniSeq LRM sample sheet to PNL35.W1, added MiniSeq checklist to PNL35-W3, and added LRM3 sample sheet to PNL35.W2. Updated PNL35.JA1 to include use of IPB.

14.5. **06/27/2024:** Updated catalog numbers and added notes for the index kits to be discontinued. Removed Microseal A from the procedure because Microseal B prevents sample loss due to evaporation better during thermocycler steps. Added *Cronobacter* and *Yersinia* with their minimum coverage requirements to Table 1. Removed all references to SPB. Added a strong recommendation to use UD indexes instead of CD indexes with this procedure. Adjusted Table 5 (compensation for sample loss during size selection) so that sample loss was only acceptable to up to 20 µl. If more sample was lost, then the library needs to be repped. Added fragment analysis as an optional quality control step. Workbooks PNL35.W1 and PNL35.W4 updated with the new Illumina DNA-RNA UD v3 indexes and current version of Generate FASTQ Modules. Sample sheets for IEM and LRM versions earlier than 3 were removed from PNL35.W1 and the NextSeq Sample Import Template was added as a tab together with instructions for its use in the appendix PNL35-3. The Metrics tab in PNL35.W1 was updated to remove calculation of coverage using versions of LRM2 or earlier. Also added highlighting samples with failing coverage in Metrics tab. Workbook PNL35.W2 was retired. Removed all references to BioNumerics.

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 16 of 21**

**15. APPROVAL SIGNATURES:**

Approved By: \_\_\_\_\_ Date: \_\_\_\_\_  
QA/QC Personnel

Approved By:   N/A   \_\_\_\_\_ Date:   NA   \_\_\_\_\_  
PulseNet Response and Outbreak Management Team Lead

Approved By: \_\_\_\_\_ Date: \_\_\_\_\_  
PulseNet WGS Technical Lead

Approved By: \_\_\_\_\_ Date: \_\_\_\_\_  
PulseNet Reference Outbreak Surveillance Team Lead

**16. APPENDICES:**

**Appendix PNL35-1. Combining PulseNet Organisms and Norovirus Amplicons on the Same Sequencing Run**

**NOTE:** *There is no need to decrease the DNA load of the bacterial DNA. The viral amplicon libraries can be spiked into PulseNet sequencing runs that are considered “full runs” for the sequencing kit chosen.*

1. Prepare the amplicons (maximum of 48 samples) as described in VGB.NV.METHOD.002 RT-PCR Protocol for Amplification of Full-Length Norovirus GI and GII Genomes.
2. Prepare the sequencing libraries from the sample amplicons as described in PNL35.
3. Pool equal volumes of libraries and dilute to 4 nM. Follow the calculations in Table 9 for diluting the pool to 4 nM.

**NOTE:** *The calculations are automated in the workbook PNL35.W1.*


	<b>Example Calculation</b>	<b>Equations</b>
<b>Pool concentration (Qubit value, ng/μl):</b>	10.6	
<b>Molarity (nM) = (conc/(660g/mol x 400 bp) x 10<sup>6</sup>):</b>	40.15	$=(\text{PoolConcentration}/(660*400))*10^6$
<b>Volume of pool to dilute (μl) = (nM)(x)=(4 nM)(50 μl):</b>	4.98	$=4*50/\text{Molarity}$
<b>Volume of RSB to dilute (μl) = 50 - volume of pool:</b>	45.02	$=50-\text{Volume of pool}$

**Table 9.** Calculations to dilute the Norovirus pool to 4 nM

4. Spike 2 μl of 4 nM amplicon library pool into the 4 nM bacterial library pool.
 

**NOTE:** *The viral amplicon spike-in does not significantly affect the bacterial library pool molarity.*
5. Proceed to sequencing on MiSeq following the instructions in PNL38.
  - 5.1. For Norovirus samples, a prefix “noro\_” needs to be added in front of the state keys to easily differentiate Norovirus samples from PulseNet samples.
  - 5.2. For index selection, make sure that all index pairs (bacterial and Norovirus) are unique to the run and have not been used in the previous run on this instrument.

**Workflow for Combined PulseNet/Norovirus Sequencing Runs**



**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 18 of 21**

**Appendix PNL35-2. Combining PulseNet Organisms and Cyclospora Amplicons on the Same Sequencing Run**

**NOTE:** *There is no need to decrease the DNA load of the bacterial DNA. The Cyclospora amplicon libraries can be spiked into PulseNet sequencing runs that are considered “full runs” for the sequencing kit chosen.*

1. Extract the Cyclospora DNA, prepare the amplicons (maximum of 24 samples) and DNA sequencing libraries as described in AMD.DR.C.001 State Cyclospora Outbreak Genotyping SOP.
2. Pool equal volumes of libraries, measure the concentration and size distribution and dilute to 4 nM. Follow the calculations in Table 10 for diluting the pool to 4 nM.

**NOTE1:** *The calculations are automated in the workbook PNL35.W1*


**NOTE2:** *The pool fragment size distribution typically varies between 300-350 bp. It is recommended determine the fragment size for every run.*

	Example Calculation	Equations
<b>Pool concentration (Qubit value, ng/μl):</b>	6.72	
<b>Molarity (nM) = (conc/(660g/mol x avg library size in bp) x 10<sup>6</sup>):</b>	25.45	$=(\text{PoolConcentration}/(660*400))*10^6$
<b>Volume of pool to dilute (μl) = (nM)(x)=(4 nM)(50 μl):</b>	7.86	$=4*50/\text{Molarity}$
<b>Volume of RSB to dilute (μl) = 50 - volume of pool:</b>	42.14	$=50-\text{Volume of pool}$

**Table 10.** Calculations to dilute the Cyclospora pool to 4 nM


3. Spike 20 μl of 4 nM Cyclospora library pool into the 4 nM bacterial library pool.
4. Proceed to sequencing on MiSeq following the instructions in PNL38.
  - 4.1. For Cyclospora sample naming in the Library Prep tab of the workbook, follow the standardized format outlined in the SOP AMD.DR.C.001, step 4.9.1.
  - 4.2. For index selection, make sure that all index pairs (bacterial and Cyclospora) are unique to the run and have not been used in the previous run on this instrument.

**Workflow for Combined PulseNet/Cyclospora Sequencing Runs**



**Appendix PNL35-3: Setting up a NextSeq Sample Import Template to Generate a NextSeq SampleSheet**

1. Ensure the data in the **Sample ID\***, **Well Position\***, and **Project** (if applicable) columns are formatted correctly on the NextSeq\_Import\_Sample\_Template tab of PNL35.W1. Data from the Library Prep tab should auto-populate these columns:
  - a. **Sample ID\*** should be in the following format: *Sample1-LabID-sequencer ID-YYMMDD*
  - b. **Well Position\*** should contain the Index plate-well location (i.e., A-A01)
  - c. **Project** data is pulled from the Project ID column on the Library Prep Tab and is optional. If this is left blank the Project in BaseSpace will default to the Run ID.
2. Save this tab (NextSeq\_Import\_Sample\_Template) as a .csv file to import into BaseSpace.
3. Log into BaseSpace, and select “**New Run**” and “**Run Planning**” from the drop-down menu as shown in Fig. 2:



**Figure 2: Selecting New Run>Run Planning to import Sample Template**

4. After entering **Run Name**, **Instrument Platform**, and **Secondary Analysis** location (Local or BaseSpace) as shown in Fig. 3, click “**Next**”.

The image shows the 'Create a Run' form in BaseSpace. The form is titled 'Create a Run' and has a 'Run Settings' section. The fields are: 'Run Name\*' with the value 'VL403-24-005' and a note that it can contain 255 alphanumeric characters, dashes, underscores, periods, and spaces; 'Run Description' with an empty text box and a note that it can contain 255 characters except square brackets, asterisks, and commas; 'Instrument Platform\*' with a dropdown menu set to 'NextSeq 1000/2000'; 'Secondary Analysis\*' with radio buttons for 'BaseSpace' (selected) and 'Local', and a note that storage and compute iCredit charges may apply; and 'Library Tube ID' with an empty text box and a note that it can contain 255 characters max.

**Figure 3: Designation of Run Settings**

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**


**Ver. No. 05**

**Effective Date:**

**Page 20 of 21**

5. Configuration (Fig. 4):

- a. Configure the application for DRAGEN BCL Convert. **Make sure this version matches what is on your instrument.**
- b. Designate the Library Prep Kit and the Index Adapter Kit used. **NOTE: Do not select individual Index Sets (i.e., Illumina DNA-RNA UD Indexes Set A Tagmentation) as this workbook generates the Well Position entry that incorporates the Index plate ID into the Well Position (A-A01). This designation is compatible with the Combined Index Adapter Kit listed in BaseSpace as “Illumina (or IDT-ILMN) DNA-RNA Indexes Set A B C D Tagmentation only. If you select the individual Index Kits for the Index Adapter Kit (i.e. Set A), this will change the format of the Import Samples Template and will give an error indicating invalid well position.**
- c. Once complete, click “Next”.




The screenshot shows the 'Create a Run' interface in BaseSpace. At the top, it displays 'BaseSpace VL403-24-005'. Below this is a 'Configuration' section with four fields:

- Application\* ?**: A dropdown menu with 'Illumina DRAGEN BCL Convert - 4.2.7' selected.
- Description**: A text input field with a '140 characters max' limit.
- Library Prep Kit\* ?**: A dropdown menu with 'Illumina DNA Prep' selected. Below the dropdown is a description: 'A fast, integrated workflow for a wide range of applications, from human whole-genome sequencing to amplicons, plasmids, and microbial species.'
- Index Adapter Kit\* ?**: A dropdown menu with 'Illumina DNA-RNA UD Indexes Set A B C D Tagmentation' selected.

**Figure 4:** Configuration of run application and kits

6. Ensure that the information for **Index Reads**, **Read Lengths**, and **Override Cycles** are correct and Select “**Import Samples**” and “**CSV**” to import the NextSeq\_Import\_Samples\_Template .csv file that was created from the NextSeq\_Import\_Samples\_Template tab of PNL35.W1 (Fig. 5).



The screenshot shows the 'Import Samples' dialog box. At the top, there are two buttons: 'Download template' and 'Import Samples'. Below the buttons is a table with the following structure:

Sample ID*	Index Reads	Read Lengths
<input type="checkbox"/>	CSV	17 index*
<input type="checkbox"/>	SAMPLE SHEET	

Below the table, there are two 'Select' buttons with dropdown arrows.

**Figure 5:** Importing NextSeq\_Import\_Samples\_Template .csv file

**PULSENET STANDARD OPERATING PROCEDURE FOR DNA LIBRARY PREPARATION USING THE ILLUMINA® DNA PREP KIT**

**Doc. No. PNL35**

**Ver. No. 05**

**Effective Date:**

**Page 21 of 21**

7. The required fields should be populated as shown below in Fig. 6.

		Sample ID*	Well Position*	I7 index*	Index 1*	I5 index*	Index 2*	Project
II	1	2013L-5214	A-C01	UDP0003	CGTCTCATAT	UDP0003	TATAGTAGCT	
II	2	2013L-5351	A-E01	UDP0005	GACGAGATTA	UDP0005	ACATTATCCT	
II	3	2013L-5356	A-F06	UDP0046	AGATCCATTA	UDP0046	ATCTCTACCA	
II	4	2013L-5357	A-E07	UDP0053	GGAATTGTAA	UDP0053	AGCACATCCT	

**Figure 6:** Confirmation of Samples/Indexes imported

8. Under **Analysis Setting**, leave **AdapterRead1 and 2** and **BarcodeMismatchesIndex1 and 2** empty as default settings (Fig. 7) and click “Next”.

**Analysis Setting**

AdapterRead1 ?

AdapterRead2 ?

BarcodeMismatchesIndex1 ?

BarcodeMismatchesIndex2 ?

Back

**Figure 7:** Analysis Settings

9. Review the run information created as shown in Fig. 8 and select **Save as Draft** or **Save as Planned**. This will create a SampleSheet in BaseSpace. You can now log into BaseSpace on the instrument, select your configured run, and load the instrument for sequencing.

**Create a Run**

**Run Review**

Run Name: test [Edit](#)

Instrument Platform: NextSeq 1000/2000

Read Lengths	Read 1	Index 1	Index 2	Read 2
	151	10	10	151

Secondary Analysis: BaseSpace

Configuration: Illumina DRAGEN BCL Convert - 4.2.7 [Edit](#)

Library Prep Kit: Illumina DNA Prep

Index Adapter Kit: Illumina DNA-RNA UD Indexes Set A B C D Tagmentation

Override Cycles: Y151;I10;I10;Y151

Samples: 66 samples

**Figure 8:** Review Run Information