

# ICLN's 2014 Full Scale Radiological Laboratory Exercise

Robert L. Jones, PhD and David Saunders, PhD

2015 APHL Meeting

## Background

After a major RDD, how will the federal, state and private radiation laboratories coordinate their analytical processes, implement surge capacity, request lab capacity from other lab networks and have a common reporting platform to have a better and more efficient and effective overall incident response to support critical decision making?

# ICLN Full Scale Radiological Lab Exercise

## ➤ Radiological Laboratory Full Scale Exercise:

- Assess Analytical Throughput, Laboratory Quality Objectives, Sample Tracking, Reporting (to networks and ICLN Portal) on real spiked samples in various matrices.
- Assess the ability to request and receive surge capacity samples
- Extensive use and testing of the **ICLN Portal** (SITREPs, requests for surge capacity, data uploading and testing of the **Minimum Data Elements format**)

# CDC'S ICLN Exercise 2014 GOALS

Test the CDC's ability, for 200 spot urine samples, to:

- receive, login, and process
- analyze (screen, identify and quantify)
- perform analytical quality control
- report analytical results
- use the new Division database to:
  - order tests
  - manage sample test and result information
  - produce QC charts/validating analytical run QC
  - produce standard CLIA compliant reports
  - generate ICLN MDE format output files
- upload MDE data files to the ICLN Data Exchange Portal

# Results - Gross $\alpha/\beta$

Spike values: Pu-239 low =  $9.2 \times 10^{-3}$  Bq/L and high =  $4.6 \times 10^{-2}$  Bq/L  
Sr-90 low = 200 Bq/L and high = 1000 Bq/L

Gross $\alpha$ (CPS/L)	Mean	SD	
Low Pu-239 $\alpha$	-0.14	1.28	n=15
High Pu-239 $\alpha$	0.16	1.24	n=15
Low Sr-90 crosstalk $\beta$ to $\alpha$	10	14	n=15
High Sr-90 crosstalk $\beta$ to $\alpha$	41	51	n=15
Blank pool $\alpha$	0.19	1.13	n=30
No Sr-90/Pu-239 $\alpha$	0.50	1.47	n=140
Gross $\beta$ (CPS/L)			
Low Sr-90 $\beta$	443	34	n=15
High Sr-90 $\beta$	2097	53	n=15
Blank pool $\beta$	37.7	7.2	n=30
No Sr-90/Pu-239 $\beta$	36.3	19.2	n=140

# Results – Pu-239 and Sr-90

Spike values: Pu-239 low =  $9.2 \times 10^{-3}$  Bq/L and high =  $4.6 \times 10^{-2}$  Bq/L

Sr-90 low = 200 Bq/L and high = 1000 Bq/L

## Pu-239, Bq/L

	<u>Mean</u>	<u>SD</u>	
Low Spike	9.04E-03	4.16E-04	n=15
High Spike	4.51E-02	1.07E-03	n=15
No Pu Spike	4.46E-06	3.71E-04	n=170
Blank pool	4.84E-05	3.62E-04	n=30

## Sr-90, Bq/L

	<u>Mean</u>	<u>SD</u>	
Low Spike	199.5	15.4	n=15
High Spike	979.9	39.9	n=15

# Issues/Lessons Identified


- Spikes thawed in Rad Team Sink water bath – cups bulged and tipped.
- Pu-239 analyte was not checked in the database analysis panel screen.
- Aliquot ID for “DLS Creatinine” was assigned during sample setup/login, but the ID for “Rad Creatinine” should have been assigned.
- “Sample login mismatch error” occurred due to use of a sample login spreadsheet that our database required be completely filled in to work. Switched to use of the “generic long form” .
- Autosampler for the Pu-239 analysis stalled.
- Pu-239 autosampler got worse and worse. Replaced it after the 5th batch.

# Issues/Lessons Identified (continued)

- Data for the 6th batch of Pu-239 samples included the wrong sample IDs. Analyst corrected this before uploading to the database.
- Gross Beta Screen and Sr-90 analysis results were processed by hand entry of raw output to pre-configured calculation spreadsheets, and there were problems uploading the data.
- Gross Beta and Sr-90 final reports and folders were initially in the wrong format.
- Creatinine QC drifted and failed at the end of an early batch. Analyst decided to run groups of 50 instead of 100.




# Pu-239 Timeline




**ABOUT 48 HOURS PER 200 SAMPLES**

# Gross Alpha/Beta Screen Timeline



**ABOUT 24 HOURS FOR 200 SAMPLES**

# Sr-90 Timeline



- Available data was not in the database, but could have been manually transferred to a report. It was imported to the database much later.
- If we had not used Gross Alpha/Beta Screening, all Sr data for complete Sr analysis of 200 samples would have been available at  $T_0 + 79$  hours (if analyzed serially).
- Cost of columns for 200 analyses would be \$4400, vs \$660 for 30.

**ABOUT 26 HOURS FOR 200 SAMPLES**

# Exercise Timeline “Take-aways”

- Manual aliquoting took a significant amount of time (7.5 to 12.5 hours) with just two people.
- We can Gross Alpha/Beta screen 200 urine samples in 24 hours (one instrument / one analyst).
- We can identify and quantify Pu-239 in 200 urine samples in 48 hours (one instrument / one analyst).
- In conjunction with Gross Alpha/Beta screening, we can identify and quantify Sr-90 in 200 urine samples (at a 15% positive Gross Alpha/Beta rate) in 26 hours.
- Gross screening prior to Sr-90 analysis can save a significant number of personnel hours and significantly reduce supply costs.

# True Emergency Timeline Modifiers

- In a real emergency we will have access to at least 40 sample logistics staff, rather than just 2, and would use an automated aliquoting system.
- Additional analysts and 4 more high capacity instruments should allow 400-500 Gross Alpha/Beta analyses in 24 hours.
- Additional analysts and 2 more Thermo Element SF-ICP-MS instruments should allow 300 Pu-239 analyses in 24 hours.
- Additional analysts, 6 more SPE processing stations, and 4 more high capacity beta analysis (LSC and GFPC) instruments should allow Sr analyses of 400-500 Gross Alpha/Beta screened samples in 24 hours.

# Summary (all)

- The Exercise has been considered a success by all the Networks since all objectives were completed.
- This is the first ICLN radiation laboratory full scale exercise evaluating surge capacity issues.
- This is the first ICLN radiation laboratory full scale exercise evaluating multiple phases of a response.
- This provided valuable Lessons Identified/Learned through the exercise process.
- This will allow for the evaluation of several ICLN Portal improvements (SitReps, data reporting, downloading, etc.) that have been recently implemented.

# Summary (CDC)

- CDC realized most of its stated goals for participating in this “Confidence Building Competency Test” (CBCT).
- CDC’s methods were effective in identifying and quantifying the radionuclides of interest.
- CDC determined that it can perform these analyses on 200 samples in 1 to 2 days.
- A number of issues were identified that can be used to improve performance in these and other methods and future exercises or responses.

# Acknowledgements

- FDA Lab Network
- EPA Lab Network
- DOE Lab Network
- Full Scale Lab Exercise Team Members
- DHS S&T ICLN
- DHS Financial Support



Questions?

# Contact

Robert L. Jones, PhD

Centers for Disease Control and Prevention

4770 Buford Hwy

Mailstop F-50

Atlanta, GA 30341-3724

RLJones@cdc.gov

“The findings and conclusions in this study are those of the authors and do not necessarily represent the views of the U.S. Department of Health and Human Services, or the U.S. Centers for Disease Control and Prevention.

Use of trade names and commercial sources is for identification only and does not constitute endorsement by the U.S. Department of Health and Human Services, or the U.S. Centers for Disease Control and Prevention.”